



Europäisches Patentamt
European Patent Office
Office européen des brevets

Publication number:

2
**0 289 896
A1**

2

EUROPEAN PATENT APPLICATION

21 Application number: 88106566.8

51 Int. Cl. 4: **C12N 5/00 , A61K 35/14**

22 Date of filing: 23.04.88

30 Priority: 27.04.87 US 42998

43 Date of publication of application:
09.11.88 Bulletin 88/45

54 Designated Contracting States:
AT BE CH DE ES FR GB GR IT LI LU NL SE

7 Applicant: **E.I. DU PONT DE NEMOURS AND
COMPANY**
1007 Market Street
Wilmington Delaware 19898(US)

Applicant: **HAEMONETICS CORPORATION**
400 Wood Road
Braintree Massachusetts 02184(US)

72 Inventor: **Dunn, George Franks, Jr.**
306 West Mantua Avenue
Wenonah New Jersey 08090(US)
Inventor: **Irr, Joseph David**
386 Briar Lane
Newark Delaware 19711(US)
Inventor: **Halpern, Lise Nadine**
17 Old Winchester Street
Newton Highlands Massachusetts 02161(US)

74 Representative: **von Kreisler, Alek,**
Dipl.-Chem. et al
Deichmannhaus am Hauptbahnhof
D-5000 Köln 1(DE)

54 A simplified method for the preparation of human lymphokine activated killer cells.

57 A lymphocyte-containing white blood cell fraction obtained by standard leukapheresis, elutriation leukapheresis or standard centrifugation can be used for production of lymphokine activated killer cells by incubation with IL-2. Removal of red blood cells and granulocytes by centrifugation on ficoll is not necessary.

EP 0 289 896 A1

A SIMPLIFIED METHOD FOR THE PREPARATION OF HUMAN LYMPHOKINE ACTIVATED KILLER CELLS

FIELD

5 This Invention pertains to adoptive immunotherapy, more particularly to the in vitro generation of human lymphokine activated killer cells for use in such therapy.

BACKGROUND

10

Adoptive immunotherapy has recently produced encouraging clinical results against some forms of cancer. See articles in the Wall Street Journal, April 9, 1987, and Time Magazine, April 20, 1987. The therapy involves removing peripheral blood from a patient, removing red blood cells (RBC's) from the blood to produce a lymphocyte-containing white blood cell (WBC) fraction, incubating the blood fraction in culture medium with interleukin-2 (IL-2) to produce activated, tumor-destroying lymphocytes called LAK cells, and injecting the LAK cells and additional IL-2 into the patient. In some cases IL-2 is injected into the patient before removal of the blood in order to stimulate production of lymphocytes.

One objection to adoptive immunotherapy is that it is very expensive. One reason it is expensive is that the current procedure for producing LAK cells is labor-intensive and time consuming. This procedure is described in Muul et al., "Large scale production of human lymphokine activated killer cells for use in adoptive immunotherapy," Journal of Immunological Methods, 88:265-275 (1986). As described in Muul et al., in order to generate enough LAK cells for a single treatment about 2×10^8 lymphocytes were obtained by 10 successive leukaphereses of peripheral blood. In each leukapheresis, about 10-12 liters of whole blood were processed in an automated cell separator over a 4-hour period to produce a 400-500 ml leucocyte fraction. This fraction was diluted with 2 parts of a salt solution, then poured into 50 ml conical centrifuge tubes (40 ml/tube, approx. 30-40 tubes) and underlayered with 10 ml Ficoll-Hypaque solution. The contents were centrifuged, causing separation into a platelet-rich supernatant layer, a lymphocyte-rich layer, a Ficoll-Hypaque layer, a granulocyte layer and RBC layer. The supernatant was removed from each tube and discarded. The lymphocyte-rich fraction floating on the Ficoll-Hypaque was removed from each tube; these fractions were pooled and washed three times by suspension in salt solution and centrifugation. Since these steps must be repeated for each leukapheresis, 300-400 tubes must be handled for a single treatment.

Haemonetics Corporation of Braintree, Massachusetts, markets an automated blood cell separator known as the Haemonetics V-50, which utilizes a 2-port conically-shaped centrifuge bowl similar to the bowl described in U.S. Patent 3,145,713. The V-50 can be operated according to a standard leukapheresis protocol or according to a Surge® lymphocytopheresis protocol. The latter procedure, as described in U.S. Patents, 4,464,167 and 4,416,654, involves intermittent elutriation with previously-separated plasma, and is capable of providing more precise fractions of platelets, WBC's and RBC's than can be achieved with standard leukapheresis; it is referred to hereinafter elutriation leukapheresis.

For LAK cell processing, Haemonetics recommends use of the V-50 to separate a Buffy coat composed mostly of platelets and WBC's, followed by a secondary separation using Ficoll-Hypaque to provide a density gradient in the same centrifuge bowl for isolation of mononuclear cells (lymphocytes and monocytes) from the Buffy coat. Although this procedure is much less time-consuming and labor-intensive than the standard ficoll centrifugation described in Muul et al., it would be desirable to eliminate the ficoll separation step because it adds to the cost and can cause a reduced yield of lymphocytes. However, up to now it has been considered essential by those skilled in the art to conduct a ficoll separation in order to obtain a lymphocyte fraction sufficiently free of RBC's and granulocytes to be useful for production of LAK cells. It was assumed that RBC's and granulocytes would unduly interfere with the activation of the lymphocytes.

50

SUMMARY OF THE INVENTION

We have discovered that the step of ficoll density gradient centrifugation can be eliminated without unduly interfering with lymphocyte activation. Thus, our invention is an improvement in the method of producing LAK cells *in vitro* which comprises removing RBC's from whole blood to produce a lymphocyte-containing WBC-rich fraction and incubating the WBC-rich fraction in culture medium with IL-2 to activate the lymphocytes. The improvement comprises using the lymphocyte-containing WBC-rich fraction without intermediate separation of a lymphocyte and monocyte layer on a ficoll gradient.

DESCRIPTION OF THE PREFERRED EMBODIMENTS

In the improved method for this invention, the RBC's can be removed in various ways. These include standard leukapheresis, elutriation leukapheresis, and centrifugation without use of ficoll. Ficoll is a synthetic water-soluble polysaccharide that has a weight average molecular weight of about 400,000 and that is widely used for the preparation of density gradients. It is available as such and in admixture with other substances under registered trademarks such as Ficoll-Paque, Ficoll-Hypaque and Ficoll-Isopaque. Leukapheresis refers to a process in which peripheral blood is withdrawn from a patient or donor, a WBC-rich fraction is separated out, and other blood fractions (plasma, platelets and RBC's) are returned to the source. Standard centrifugation is used to separate blood from donors into plasma, WBC-rich and RBC fractions which are stored for later use. (The term "Buffy coat" as used hereinafter refers specifically to the WBC-rich fraction obtained by standard centrifugation, although the term is also used in the art to refer to a platelet-rich, WBC-rich fraction obtained by leukapheresis.)

The various methods of removing RBC's produce WBC-rich fractions with varying amounts of residual RBC's and varying differentials. (The term "differential" or "diff" refers to the number percent of lymphocytes, monocytes and granulocytes based on the total number of those three cell types in a WBC-rich fraction.) Compositions of the various fractions will also vary depending upon the source. For example, a patient who has been primed with IL-2 may have a very high lymphocyte count. Typical ranges for the WBC-rich fractions obtained by various methods are compared with typical ranges for whole blood in the following table.

	No RBC	No WBC	Differential		
			L	M	G
Standard Leukapheresis per 240 ml pack Vol. % RBC 10-20	10^{11} to 5×10^{11}	2×10^9 to 10^{11}	60-80	5-25	5-25
Elutriation Leukapheresis per 400 ml pack Vol. % RBC 1-6	2×10^{10} to 10^{11}	2×10^9 to 10^{11}	80-85	10-20	1-5
Buffy Coat per 40 ml pack Vol. % RBC 40-50	10^{11} to 3×10^{11}	10^8 to 2×10^{10}	20-50	10-30	20-50
Normal Whole Blood per 450 ml unit	$1.6-2.4 \times 10^{12}$	$2.3-4.6 \times 10^9$	25-40	4-10	50-65

From the above table, it can be seen that lymphocyte-containing WBC-rich fractions usable in this invention can have RBC/WBC ratios from about 0.2 to about 300 and granulocyte contents from about 1% to about 50%. As a practical matter Buffy coats would principally be used for screening to determine whether a patient is capable of developing LAK cells. For generating LAK cells for use in adoptive

immunotherapy the leukapheresis products having RBC vol. % of about 1-20% and RBC/WBC ratio of about 0.2-250 would be preferred.

At the present time, it is preferred to use the elutriation leukapheresis product because it is more nearly like the ficoll-separated products in both RBC and granulocyte content, and therefore would probably be more readily accepted by workers in the art. In addition, it appears that elutriation leukapheresis products can be cultured at a somewhat higher cell density (e.g., 1×10^7 /ml or higher) than can standard leukapheresis products on a routine basis. From the above table, it can be seen that elutriation leukapheresis products typically have a RBC/WBC ratio of about 0.2 to about 50, a RBC vol. percent of about 1-6 and a granulocyte content of 1-5. More typical ranges are RBC/WBC of about 0.5-25 and RBC vol. % of about 2-4.

Standard leukapheresis can be performed using instruments available from various manufacturers, including Haemonetics, Fenwall, and Cobe and following the manufacturers' instructions. The only instrument now available for performing elutriation leukapheresis is the Haemonetics V-50. Following the teaching of U.S. Patents 4,464,167 and 4,416,654 or the instructions provided by Haemonetics, the V-50 can be used to provide a WBC-rich fraction having low RBC and granulocyte content.

Monocyte content of the WBC-rich fraction can be reduced below the figures shown in the table by treatment of the leukapheresis product with an L-amino acid lower alkyl ester or hydrogen chloride salt thereof, e.g., methyl, ethyl, propyl, isopropyl, butyl, isobutyl, or t-butyl ester of phenylamine, glutamic acid, glutamine or tyrosine. Phenyl alanine methyl ester is preferred. Further details are given in copending U.S. application Serial No. 868,697, filed May 30, 1986, and in the examples below.

Activation of the lymphocytes by incubation with IL-2 is accomplished in this invention in the same manner as in the prior art. Containers such as conventional flasks and roller bottles can be used, but the preferred containers are 0.2-5 liter tissue culture bags made from flexible copolymeric film materials as disclosed in copending application Serial No. 008,273, filed January 29, 1987. Most preferred is a bag made of a copolymer of 97 mol % ethylene and 1-octene. Any suitable culture medium can be used, but the preferred culture medium is RPMI 1640, which is described in "Culture of Animal Cells", Freshney, 72-73, Alan R. Liss, Inc., NY, supplemented with serum. Initial cell concentrations of up to about 1×10^7 cells/ml can be used with an elutriation leukapheresis product and up to about 1×10^7 cells/ml with a standard leukapheresis product. A concentration of at least 1×10^6 cells/ml should be used for reason of economy. Preferred ranges would be 5×10^5 to 10^7 cells/ml for elutriation leukapheresis products and $1-5 \times 10^5$ cells/ml for standard leukapheresis products. The cells are incubated with IL-2 for about 2-7 days, preferably about 3-5 days, at a temperature of about 35-39°C, preferably 37°C.

"Interleukin-2" (IL-2) as used herein means human IL-2. It includes natural and recombinant IL-2 (rIL-2) and biologically functional equivalents thereof, such as the rIL-2 mutants disclosed in U.S. Patent 4,518,584. Preferably, the IL-2 is a rIL-2 composition consisting essentially of water, rIL-2 and, optionally, a polyol as described in assignee's copending application Serial No. 825,133, filed on January 31, 1986. Preferably, the IL-2 concentration in the culture medium is in the range of about 5×10^2 to about 5×10^4 pM, most preferably 1000 to 2000 pM.

The LAK cells prepared by the process of the invention can be suspended in a pharmaceutically acceptable carrier, such as saline, saline containing 5% normal human serum albumin, or Hank's balanced salt solution, to provide a composition which can be infused into a patient afflicted with a tumor. The patient is concurrently treated with rIL-2 as further described by Rosenberg et al., The New England Journal of Medicine 313, 1485-1492 (1985). In that modality, the patient's blood is withdrawn, subjected to leukapheresis and harvested cells are immediately cultured for 3 days to generate LAK cells. The LAK cells are then infused into the patient. Typically, about 3×10^{10} to 14×10^{10} LAK cells are infused in 4-9 doses. Interleukin-2 is administered every eight hours at doses such as 10,000, 30,000 or 100,000 units per kilogram of weight. The treatment consists of a two-week regime of leukapheresis and reinfusion and generally repetition starting the third week. Recombinant IL-2 can be included in the LAK cell composition.

Cytotoxicity (LAK) Assay

In the following examples, unless otherwise stated, a 4 hour ^{51}Cr release assay was used to measure cytotoxicity of LAK cells for tumor cells (LAK activity). Tumor cells at a concentration of about 2×10^5 to 10×10^5 per ml were incubated with 100 μCi of $\text{Na}_2^{51}\text{CrO}_4$ in 0.4 mL of Tris-phosphate buffered saline for 1 hour at 37°C. The cells were washed 3 times with RPMI 1640 containing 5% or 10% fetal calf serum (FCS) and resuspended to 10^5 cells/mL in RPMI-20% FCS or RPMI-10% FCS. The effector cells (LAK

cells) were suspended to various concentrations of 0.1 mL was added in to wells round bottom microliter plates. The ^{51}Cr labelled target cells (0.1 mL) were added to all wells. After 4 hours of incubation at 37°C, the plates were centrifuged and 0.1 mL of resulting supernatant was removed from each well and counted in a gamma counter. Percent cytolysis is calculated from the following formula:

$$\% \text{ cytolysis} = \frac{\text{experimental cpm} - \text{spontaneous cpm}}{\text{total cpm} - \text{spontaneous cpm}} \times 100$$

Each variable was tested in triplicate and the resulting data are expressed as % cytolysis. This cytotoxicity test is further described in "Selected Methods in Cellular Immunology," Mishell and Shiigi, eds., 124-137, W. H. Freeman and Co., San Francisco (1980).

In other experiments, the results of the assays are presented as "Lytic Units" (LU or LU30) which are the number of target cells per 100 effector cells when 30% of the target cells are killed when LAK cells and target cells are incubated together for 4 hours at 37°C. The calculation of LU is based upon the method of Pross et al., Journal of Immunological Methods 68, 235-249 (1984). The greater the number of LU, the greater the potency of the LAK cell preparation.

All patents, patent applications and other printed publications cited in this application are incorporated herein by reference, especially the disclosure of U.S. Patents 4,464,167 and 4,416,654 relating to the production of a WBC-rich fraction by elutriation leukapheresis using previously separated plasma as elutriant.

EXAMPLE 1

Purpose:

- 1) To study the LAK activity of cells obtained from a Haemonetics V50 using the elutriation technique to obtain white blood cells.
- 2) To study the effects of phenyl alanine methyl ester (α Ala) treatment and Ficoll treatment on the LAK activity of cells obtained from the Haemonetics V50 elutriation technique.

Cells:

Human lymphocytes (obtained from Haemonetics Corporation using V50 elutriation technique). Raji cells.

Materials:

- 1) Cell culture medium (CCM) = RPMI 1640 with 10% FBS, L-glutamine and Gentamicin
- 2) Phosphate buffered saline (PBS) 1x without Ca^{++} and Mg^{++}
- 3) Ficoll Hypaque (Ficoll)
- 4) α Ala
- 5) Unopette® for WBC count
- 6) Ethylene butene copolymer Bag for cell culture
- 7) T25 tissue culture flasks
- 8) 1% NP 40
- 9) 2x TD buffer
- 10) ^{51}Cr (as sodium chromate)
- 11) recombinant Interleukin-2, 10 units/ml in 0.5M glucose (IL-2)
- 12) 96 well μ bottom tissue culture plate

- 13) SCS Harvesting System (Skatron)
- 14) Beckman Gamma 4000 Counter
- 15) Trypan Blue

5

Procedure:A) Preparation of Cells

10

- 1) A total of 250 ml of a white blood cell fraction was collected from a Haemonetics V50 machine using the elutriation procedure as described in U.S. Patents 4,416,654 and 4,464,167.
- 2) A WBC count was performed using a Unopette®. The fraction contained 1.36×10^7 WBC/ml and was estimated to contain approximately 3 vol. % RBC.
- 3) The cells were then brought to a concentration of 1×10^7 WBC/ml (total volume = 340 ml).
- 4) 40 ml of cells were put directly into culture (as described below).
- 5) The remaining 300 ml were treated with ϕ Ala (as described below).

20

B) ϕ Ala Treatment

- 1) Place 300 ml of cells into T150 flask.
- 2) Add 30 ml of ϕ Ala to cells
- 3) Mix well (gently).
- 4) Incubate at room temperature for 40 minutes.
- 5) After incubation, separate blood into 2 aliquots each containing 165 ml
 - a) aliquot 1 was placed into culture
 - b) aliquot 2 was underlayered with Ficoll (as described below) and then placed into culture.

30

C) Set Up Culture

35

1) Cells Straight from V50 (No Ficoll; No ϕ Ala)

- a) place 40 ml cells into a 50 ml centrifuge tube
- b) centrifuge cells 10 minutes at 1200 rpm
- c) discard supernatant
- d) resuspend cells in CCM to a total volume of 40 ml
- e) place desired amount of cells into T25 flasks
- f) add CCM to flasks to bring white cells to desired concentration
- g) add 5 μ l IL-2 to each flask (final concentration 10 units/ml)
- h) place flasks in 37°C incubator with 5% CO₂.

45

Set up 3 - T25 Flasks

50

<u>1×10^6 WBC/ml</u>	<u>5×10^6 WBC/ml</u>	<u>1×10^7 WBC/ml</u>
1 ml cells	5 ml cells	10 ml cells
9 ml media (CCM)	5 ml media (CCM)	5 μ l IL-2
5 μ l IL-2	5 μ l IL-2	

55

Aliquot 1 → Cells from V50 (øAla and No Ficoll)

- 5 a) place 165 ml of øAla treated cells into a 250 ml centrifuge tube
 b) centrifuge for 10 minutes at 1200 rpm
 c) discard supernatant
 d) resuspend cells in 50 ml CCM
 e) perform cell count using Unopette®; the WBC count was 1.9×10^7 per ml
 f) set up cultures in bags and flasks according to cell concentrations desired
 10 g) place cultures in 37°C incubator with 5% CO₂.

Set Up 2 Cultures: T25 Flask 5×10^6 cells/ml
 Bag 9×10^6 cells/ml

$$\frac{5 \times 10^6}{1.9 \times 10^7} = .260 \times 10 = 2.6 \text{ ml cells}$$

$$\left. \begin{array}{l} 7.4 \text{ ml media (CCM)} \\ 5 \mu\text{l IL-2} \end{array} \right\} \text{ in T25 flask}$$

$$\frac{9 \times 10^6}{1.9 \times 10^7} = .474 \times 100 = 47.4 \text{ ml cells}$$

$$\left. \begin{array}{l} 52.6 \text{ ml media (CCM)} \\ 50 \mu\text{l IL-2} \end{array} \right\} \text{ in Bag}$$

3) Aliquot 2 → Cells from V50(øAla and Ficoll)

- 35 a) place 40 ml of øAla treated cells into 4-50 ml centrifuge tubes
 b) underlay blood with 10 ml of Ficoll
 c) centrifuge for 30 minutes at 1900 rpm
 d) collect interface layer with a sterile pasteur pipette and place cells into a sterile 50 ml centrifuge
 tube
 40 e) bring volume in the tube up to 50 ml using PBS
 f) centrifuge for 10 minutes at 1200 rpm
 g) discard supernatant
 h) resuspend pellet in 50 ml of CCM
 i) centrifuge for 10 minutes at 1200 rpm
 45 j) resuspend in 5 ml of CCM
 k) perform cell count using trypan blue: then
 l) set up cultures in bags and flasks according to cell concentrations desired
 m) place cultures in 37°C incubator with 5% CO₂.

NOTE: No interface layer resulted after step 3; therefore, the cells were resuspended, re-underlayered with Ficoll and recentrifuged. After this, the cells in the interface were collected.

Set Up 3 Cultures: T25 Flask 5×10^6 cells/ml
 T25 Flask 10×10^6 cells/ml
 Bag 1.9×10^6 cells/ml

$$\frac{5 \times 10^6}{6.8 \times 10^7} = .074 \times 10 = .740 \text{ ml cells}$$

$$\left. \begin{array}{l} 9.260 \text{ ml media (CCM)} \\ 5 \text{ } \mu\text{l IL-2} \end{array} \right\} \text{ in T25 flask}$$

$$\frac{1 \times 10^7}{6.8 \times 10^7} = .147 \times 10 = 1.47 \text{ ml cells}$$

$$\left. \begin{array}{l} 8.53 \text{ ml media (CCM)} \\ 5 \text{ } \mu\text{l IL-2} \end{array} \right\} \text{ in T25 flask}$$

$$\frac{1.9 \times 10^6}{6.8 \times 10^7} = .0279 \times 100 = 2.79 \text{ ml cells}$$

$$\left. \begin{array}{l} 97.21 \text{ ml media (CCM)} \\ 50 \text{ } \mu\text{l IL-2} \end{array} \right\} \text{ in T25 flask}$$

D) LAK Assay

The LAK assay was performed after cells were in culture for 4 days, according to the procedure given above.

NOTE: Due to the overabundance of red blood cells contained in the specimens, 3 specimens were treated with lysis buffer prior to the LAK assay.

Lysing solution:

.83 g Ammonium Chloride

200 ml distilled H_2O

- 1) resuspended cell pellet in 10 ml lysing solution
- 2) incubate for 20 minutes at room temperature
- 3) centrifuge for 10 minutes at 1200 rpm
- 4) discard supernatant
- 5) resuspend in 1 ml of CCM
- 6) perform cell count
- 7) set up E:T ratios as described in LAK assay procedure.

DATACell Counts

<u>Specimen</u>	<u>Counts</u>	<u>Cells</u>	<u>Total Cells</u>
5 Straight from V50	68	1.36×10^7	$3.4 \times 10^9 / 250 \text{ ml}$
After ϕ Ala	97	1.9×10^7	$9.7 \times 10^8 / 50 \text{ ml}$
10 After ϕ Ala and Ficoll	344	6.8×10^7	$3.4 \times 10^8 / 5 \text{ ml}$

Prepare Cells for ϕ Ala Treatment (Bring all cells to
 1×10^7 cells/ml)

$$\frac{1 \times 10^7}{1.36 \times 10^7} = .735 \times 340 = 250 \text{ ml cells}$$

90 ml CCM

20

Take off 40 ml and put into culture.
Add 30 ml of ϕ Ala to remaining cells and incubate 40 minutes at room temperature.

25

30

35

40

45

50

55

Day 4 - ⁵¹Cr Release Data

<u>Specimen</u>	<u>Viable</u>	<u>Non-Viable</u>	<u>% Viable</u>	<u>Cells/ mls</u>	<u>Mls of Cells</u>	<u>Mls of Media</u>	<u>Total Cells</u>
<u>Straight from V50</u>							
*1x10 ⁶ Flask	39	3	92%	7.8x10 ⁶	.256	.744	
*5x10 ⁶ Flask	113	44	87%	2.2x10 ⁷	.091	.909	
1x10 ⁶ Flask	136			2.7x10 ⁷	.074	.926	5.4x10 ⁷
<u>4Ala No Ficoll</u>							
*5x10 ⁶ Flask	44	2	95%	8.8x10 ⁶	.228	.772	
F9x10 ⁶ Bag	92			1.8x10 ⁷	.111	.889	3.6x10 ⁷
<u>4Ala and Ficoll</u>							
5x10 ⁶ Flask	154	21	88%	3x10 ⁷	.067	.933	
10x10 ⁶ Flask	158	30	84%	3.2x10 ⁷	.063	.937	6.4x10 ⁷
1.9x10 ⁶ Bag	48	2	92%	9.6x10 ⁶	.208	.792	
Raji	82	5	94%	1.6x10 ⁷	.240	39.760	

Straight From V50

1x10 ⁶ Flask	$\frac{2 \times 10^6}{7.8 \times 10^6} = .256$
5x10 ⁶ Flask	$\frac{2 \times 10^6}{2.2 \times 10^7} = .091$
10x10 ⁶ Flask	$\frac{2 \times 10^6}{2.7 \times 10^7} = .074$

Ala No Ficoll

$$5 \times 10^6 \text{ Flask } \frac{2 \times 10^6}{8.8 \times 10^6} = .228$$

$$9 \times 10^6 \text{ Bag } \frac{2 \times 10^6}{1.8 \times 10^7} = .111$$

Ala and Ficoll

$$1 \times 10^6 \text{ Flask } \frac{2 \times 10^6}{3 \times 10^7} = .067$$

$$10 \times 10^6 \text{ Flask } \frac{2 \times 10^6}{3.2 \times 10^7} = .063$$

$$1.9 \times 10^6 \text{ Bag } \frac{2 \times 10^6}{9.6 \times 10^6} = .208$$

$$\text{Raji} = \frac{1 \times 10^5}{1.6 \times 10^7} = .006 \times 40 = .240 \text{ ml cells} \\ 39.760 \text{ ml media}$$

* specimens which were processed with lysing solution
prior to LAK assay

F specimens counted using Unopette* method

Day 4 - ⁵¹Cr Release Data

Maximum (Total) Release = 2011	Spontaneous Release = 463	21.0%
2215	469	21.3%
<u>2361</u>	<u>573</u>	<u>26.0%</u>
Avg = 2196 cpm	Avg = 502cpm	22.8%

Cells Straight From V50 - No δ Ala and No Ficoll

Dilution	Flask 1x10 ⁶		Flask 5x10 ⁶		Flask 1x10 ⁷	
	CPM	% Cyto- lysis	CPM	% Cyto- lysis	CPM	% Cyto- lysis
20:1	1011	30.1	1102	35.4	1222	42.5
	978	28.1	1143	37.9	1437	55.2
	1068	33.4	1142	37.8	1321	48.4
	mean	<u>30.5</u>	mean	<u>37.0</u>	mean	<u>48.7</u>
10:1	821	18.9	752	14.8	1014	30.2
	736	13.8	861	21.2	1072	33.7
	758	15.1	861	21.2	1361	50.7
	mean	<u>15.9</u>	mean	<u>19.1</u>	mean	<u>38.2</u>
5:1	613	6.6	757	15.1	907	23.9
	638	8.0	719	12.8	874	22.0
	591	5.3	747	14.5	815	18.5
	mean	<u>6.6</u>	mean	<u>14.1</u>	mean	<u>21.5</u>
2.5:1	512	0.6	635	7.9	584	4.9
	548	2.7	699	11.6	629	7.5
	495	-0.4	592	5.3	1012	30.1
	mean	<u>1.0</u>	mean	<u>8.3</u>	mean	<u>14.2</u>

 δ Ala and No Ficoll

Dilution:	Flask 5x10 ⁶		Bag 9x10 ⁶	
	CPM	% Cyto- lysis	CPM	% Cyto- lysis
20:1	838	19.9	670	9.9
	716	12.7	962	27.2
	727	13.3	1158	38.7
	mean	<u>15.3</u>	mean	<u>25.3</u>
10:1	624	7.2	811	18.3
	550	2.9	723	13.1
	636	7.9	764	15.5
	mean	<u>6.0</u>	mean	<u>15.6</u>
5:1	541	2.3	719	12.8
	506	0.3	713	12.5
	654	9.0	822	18.9
	mean	<u>3.9</u>	mean	<u>14.7</u>
2.5:1	476	-1.5	709	12.2
	488	-0.8	713	12.5
	405	-5.7	727	13.3
	mean	<u>-2.7</u>	mean	<u>12.7</u>

øAla and No Ficoll

5	<u>Dilution</u>	<u>Flask 1x10⁶</u>		<u>Flask 5x10⁶</u>		<u>Flask 1x10⁷</u>	
		<u>CPM</u>	<u>% Cyto- lysis</u>	<u>CPM</u>	<u>% Cyto- lysis</u>	<u>CPM</u>	<u>% Cyto- lysis</u>
10	20:1	1030	31.2	1239	43.5	1498	58.8
		1010	30.0	1200	41.2	1227	42.8
		1069	33.5	1288	46.4	1187	40.5
		mean	<u>31.6</u>	mean	<u>43.7</u>	mean	<u>47.4</u>
15	10:1	769	15.8	982	28.4	1125	36.8
		781	16.5	999	29.4	1188	40.5
		745	14.4	939	25.8	1261	44.8
		mean	<u>15.5</u>	mean	<u>27.8</u>	mean	<u>40.7</u>
20	5:1	623	7.2	876	22.1	1127	36.9
		624	7.2	755	15.0	953	26.6
		676	10.3	759	15.2	1326	48.7
		mean	<u>8.2</u>	mean	<u>17.4</u>	mean	<u>37.4</u>
25	2.5:1	422	-4.7	627	7.4	983	28.4
		464	-2.2	556	3.2	940	25.9
		426	-4.5	563	3.7	894	23.2
		mean	<u>-3.8</u>	mean	<u>4.8</u>	mean	<u>25.8</u>

30

35

40


45

50

55

EXAMPLE 2

Protocol: The following diagram summarizes the protocol for this example.

Procedure:

A) Separation of Cells

- 1) Cells were collected via elutriation technique on Haemonetics V-50.
- 2) A cell count was performed = 1.3×10^7 cells/ml in 442 mls. 5.8×10^9 total cells.
- 3) Cell Volume was split in two for processing.

B) LAB Ficoll

- 1) 221 mls of cells were mixed with PBS and layered on Ficoll.
- 2) Centrifuged 30 min. at 2000 rpm
- 3) Cells were then washed and counted
- 256 cells viable
- 99% Viability
- 1 nonviable cell
- 5.1×10^8 cells ml
- $2 \times 10^9/40$ ml
- 4) Set up cells in Culture for LAK @ 1×10^5

Sample Code #	Conc. Wanted	Calculation	ml cells + ml Media + μ l IL-2
1	100		
2	100		
3	100		
4	100		
5	100		
6	100		
7	100		
8	100		
9	100		
10	100		
11	100		
12	100		
13	100		
14	100		
15	100		
16	100		
17	100		
18	100		
19	100		
20	100		
21	100		
22	100		
23	100		
24	100		
25	100		
26	100		
27	100		
28	100		
29	100		
30	100		
31	100		
32	100		
33	100		
34	100		
35	100		
36	100		
37	100		
38	100		
39	100		
40	100		
41	100		
42	100		
43	100		
44	100		
45	100		
46	100		
47	100		
48	100		
49	100		
50	100		
51	100		
52	100		
53	100		
54	100		
55	100		
56	100		
57	100		
58	100		
59	100		
60	100		
61	100		
62	100		
63	100		
64	100		
65	100		
66	100		
67	100		
68	100		
69	100		
70	100		
71	100		
72	100		
73	100		
74	100		
75	100		
76	100		
77	100		
78	100		
79	100		
80	100		
81	100		
82	100		
83	100		
84	100		
85	100		
86	100		
87	100		
88	100		
89	100		
90	100		
91	100		
92	100		
93	100		
94	100		
95	100		
96	100		
97	100		
98	100		
99	100		
100	100		

5 #4 1×10^6 $\frac{1 \times 10^6}{5.1 \times 10^7 \times 10} = 0.2 \text{ ml} + 9.8 \text{ ml} + 5 \mu\text{l}$

5) The remaining of these Ficoll layered cells were set up for ϕ Ala

$$a) V_1 C_1 = V_2 C_2$$
$$10 \quad (40 \text{ mls})(5.1 \times 10^7) = V_2 (1 \times 10^7)$$

#mls total = V_2 = 200 mls

$$\# \text{mls media} = V_2 - V_1 = 160$$
$$\# \text{ mls } \phi \text{Ala} = V_2/9 = 200/9 = 22.2 \text{ mls } \phi \text{Ala}$$

b) Incubate 40 min and then wash.

15 c) Perform cell count and put cells in culture for LAK

Viable cells = 229

Nonviable = 6

% Viability = 97%

Cells/ml = 4.6×10^7

20 d) Set up cells for culture

Sample	Conc.
Code #	Wanted Mls Cells + Ml media + μ l IL-2

25 #1 1×10^6 0.20 ml + 9.8 mls + 5 μ l

#2 5×10^6 10.9 ml + 89.1 ml + 50 μ l

#3 1×10^7 21.8 ml + 78.2 ml + 50 μ l

30

C) LAK directly from V-50

35

1) The second half of cells were used at this time. The amount of cells necessary to have cultures at a concentration of 1×10^6 , 5×10^6 and 1×10^7 were used and the remaining cells were diluted and treated with phenyl alanine methyl ester.

sample Code #	Conc. Wanted	ml Cells +	ml media +	ul IL-2
40				

#5 1×10^6 0.76 ml + 9.2 mls + 5 μ l

#6 5×10^6 3.9 ml + 6.1 ml + 50 μ l

45	#6	5x10	3.9 ml	+ 6.1 ml	+ 50 μ l
	#7	1x10 ⁷	7.7 ml	+ 2.3 ml	+ 50 μ l

#7A 1×10^7 7.7 ml + 2.3 ml + 50 μ l added

*10 ml sample - centrifuged - removed 5 mls plasma;

50 added 5 mls media. Diluted cells to 1×10^7 .

55

D) ϕ Ala without Ficoll

ϕ Ala cells without separating with Ficoll first

$$C_1 V_1 = C_2 V_2$$

$$(200 \text{ ml})(1.2 \times 10^7) = (x)(1 \times 10^7) = x = 260$$

mls media = 60 mls

mls ϕ Ala = 29 mls.

Incubate 40 min.

10 Wash.

Perform cell count and put up in culture

Viable 240

%Viable 93%

Cells/ml = 4.8×10^7

15 Total = $1.9 \times 10^9 / 40 \text{ ml}$

Sample Code #	Conc. Wanted	ml Cells + ml media + μ l IL-2
#8	1×10^6	0.2 ml + 9.8 mls + 5 μ l
#9	5×10^6	10.4 ml + 89.6 ml + 50 μ l
#10	1×10^7	20.8 ml + 72.2 ml + 50 μ l

25 All cultures were incubated for 3 days at 37°C and 5% CO₂. LAK - ⁵¹Cr release performed. In the tables which follow, % REL means % Release, which is the same as % Cytolysis, calculated as explained above. E:T means the ratio of effector (LAK) cells to target (tumor) cells.

Cell counts were performed on all cultures.

30

35

40

45

50

55

	Sample	Viable	Non-Viable	% Viable	Total
	A	33	4	88	6.5×10^6 /ml
	B	99	8	93	2.0×10^8 /10 ml
5	C	236	17	93	4.7×10^8 /10 ml
	D	7	0	100	1.4×10^6 /ml
	*E	22	1	96	4.4×10^6 /ml
10	F	45	5	91	1.8×10^7 /2 ml
	G	38	4	90	7.6×10^6 /ml
	H	75	6	93	1.5×10^8 /10 ml
15	J	141	15	90	2.8×10^8 /10 ml
	K	31	4	89	6.2×10^6 /ml
	# 1	34	1	97	6.8×10^6 /ml
20	2	162	8	95	3.2×10^8 /10 ml
	3	375	20	94	7.5×10^8 /10 ml
	4	41	2	95	8.1×10^6 /ml
25	5	27	2	93	5.4×10^6 /ml
	6	43	1	98	1.7×10^7 /2 ml
	7	121	17	88	4.8×10^7 /2 ml
30	7A	149	9	94	6.0×10^7 /2ml
	8	50	0	100	1×10^7 /ml
	9	113	5	96	2.3×10^8 /10 ml
35	10	424	47	90	8.5×10^8 /10 ml

*Letter E had a heavy fibrin clot after centrifugation.
 Raji tumor cells were prepared as target - Viability
 was 100% with a concentration of 3.6×10^5 /ml

#1 @ 1×10^6

P 05

19.77

0 289 896

#3 @ 1×10^7

	<u>E:T</u>	<u>CPM</u>	<u>% REL</u>
	20	603	76.34
5	20	647	83.72
	20	571	70.97
	10	546	66.78
10	10	408	43.62
	10	405	43.12
	5	351	34.06
15	5	333	31.04
	5	324	29.53
	2.5	265	19.63
	2.5	217	11.58
20	2.5	271	20.64

STD Ficoll - No 6Ala

#4 @ 1×10^6

	<u>E:T</u>	<u>CPM</u>	<u>% REL</u>
	20	565	69.97
30	20	507	60.23
	20	529	62.42
	10	483	56.21
35	10	400	42.28
	10	462	52.68
	5	276	21.48
40	5	255	17.95
	5	291	23.99
	2.5	288	23.49
45	2.5	263	19.30
	2.5	235	14.60

50

55

Elutriation - No Ficoll or ϕ Ala#5 @ 1×10^6

	<u>E:T</u>	<u>CPM</u>	<u>% REL</u>
5	20	588	73.83
	20	496	58.39
	20	466	53.36
10	10	305	26.34
	10	296	24.83
	10	339	32.05
15	5	285	22.99
	5	281	22.32
	5	228	13.42
20	2.5	211	10.57
	2.5	229	13.59
	2.5	237	14.93

25

#6 @ 5×10^6

	<u>E:T</u>	<u>CPM</u>	<u>% REL</u>
	20	759	102.52
30	20	626	80.20
	20	512	61.07
	10	543	66.28
35	10	533	64.60
	10	521	62.58
	5	373	37.75
	5	474	54.70
40	5	408	43.62
	2.5	224	12.75
	2.5	341	32.38
45	2.5	266	19.80

50

55

#7 @ 1×10^7

	<u>E:T</u>	<u>CPM</u>	<u>% REL</u>
5	20	314	27.85
	20	288	23.49
	20	312	27.52
	10	203	9.23
10	10	186	6.38
	10	206	9.73
	5	187	6.54
15	5	185	6.21
	5	177	4.87
	2.5	157	1.51
	2.5	190	7.05
20	2.5	167	3.19

#7A @ 1×10^7

	<u>E:T</u>	<u>CPM</u>	<u>% REL</u>
25	20	481	55.87
	20	308	26.85
	20	297	25.00
30	10	249	16.95
	10	281	22.32
	10	252	17.45
	5	197	8.22
35	5	144	-0.67
	5	167	3.19
	2.5	160	2.01
40	2.5	197	8.22
	2.5	164	2.68

45

50

55

Elutriation - 6Ala. No Ficoll

#8 @ 1×10^6

	<u>E:T</u>	<u>CPM</u>	<u>% REL</u>
5	20	592	74.50
	20	574	71.48
	20	463	52.85
10	10	377	38.42
	10	438	48.66
	10	425	46.48
15	5	257	18.29
	5	228	13.42
	5	431	47.48
20	2.5	207	9.90
	2.5	217	11.58
	2.5	211	10.57

25

#9 @ 5×10^6

	<u>E:T</u>	<u>CPM</u>	<u>% REL</u>
	20	687	90.44
30	20	573	71.31
	20	598	75.50
	10	636	81.88
35	10	609	77.35
	10	647	83.72
	5	255	17.95
	5	284	22.82
40	5	352	34.23
	2.5	264	19.46
	2.5	346	33.22
45	2.5	315	28.02

50

55

		#10 @ 1×10^7	
	<u>E:T</u>	<u>CPM</u>	<u>REL</u>
5	20	572	71.14
	20	586	73.49
	20	493	57.89
10	10	349	33.72
	10	288	23.49
	10	254	17.79
	5	183	5.87
15	5	204	9.40
	5	233	14.26
	2.5	205	9.56
20	2.5	214	11.07
	2.5	195	7.89

25

EXAMPLE 3

30 A standard leukapheresis product containing 225 mls human leukocytes prepared from 3600 mls whole blood collected in 550 mls anticoagulant ACD-B was obtained from Biological Speciality Corporation, Lansdale, PA. The following procedures were performed using this product.

1) Set up a Unopette® (WBC) and a differential. Differential:

90% Lymphocytes

6% Monocytes

35 4% Granulocytes

Direct = 165 cells

3.3×10^7 cells/ml

7.4×10^9 cells/225 mls

2) Set up cells in culture for LAK

40 @ 1×10^6 = 0.30 ml cells + 9.70 ml media + 5 μ l IL-2

@ 5×10^6 = 1.52 ml cells + 8.48 ml media + 5 μ l IL-2

@ 1×10^7 = 3.04 ml cells + 6.96 ml media + 5 μ l IL-2

Incubated @ 37°C, 5% CO₂ for 4 days.

45 3) Next 20 mls of cells were removed from the remaining cells and mixed with 20 mls of PBS without Ca⁺⁺ and Mg⁺⁺. 40 mls of this mixture was layered onto 40 mls of Ficoll and centrifuged for 1.2 hour. Removed mononuclear cell layer and washed these cells 3 times. Performed 90 min. monocyte adherence. Washed 2 more times and performed cell count and put in culture for LAK. This is the Standard Sample.

Standard Cell Count:

50 Viable = 155

Non-viable = 1

% Viable = 99%

Cells/ml = 3.1×10^7 ml

Total = 7.7×10^8 /25 ml.

55 Dilution for cell conc. of 1×10^5 = 0.32 ml cells + 9.68 ml media + 5 μ l (IL-2).

4) The remaining cells (200 mls) were then diluted with 460 mls CCM to bring the cell count to 1×10^7 ml, and treated with 73 mls ϕ Ala.

Incubated at R.T. for 40 min.

Cells clotted during incubation.

Removed as much unclotted suspension as possible --washed and counted cells with WBC Unopette®.

26 cells viable

1 cell nonviable

96% viability

5.2x10⁶ cells/ml

1.04x10⁹ cells/200 mls.

Put cells up in a bag at 5x10⁶ = 48 mls cells + 2 ml media + 25 µl IL-2.

Incubated at 37°C 5% CO₂ for 4 days for LAK.

10

5) Cell Counts After Incubation

		Viab	Non-Vi	% Viab	Cells/ml	Total
15	Std.	25	1	96	4.9x10 ⁶	---
	Direct @ 1x10 ⁶	27	1	96	5.3x10 ⁶	---
	Direct @ 5x10 ⁶	11	1	92	2.1x10 ⁶	1.1x10 ⁷ /5ml
20	Direct @ 1x10 ⁷	24	5	83	4.8x10 ⁶	4.8x10 ⁷ /10ml
	Ala @ 5x10 ⁶	19	0	100	3.8x10 ⁶	
	Raji	85	7	92	1.7x10 ⁷	

25

E:T ratio 20:1, 10:1, 5:1, 2.5:1

Dilutions

30

2x10⁶

35

Std.	0.41 mls cells + 0.59 mls media
Direct 1x10 ⁶	0.38 ml cells + 0.62 ml media
Direct 5x10 ⁶	0.95 ml cells + 0.05 ml media
Direct 11x10 ⁷	0.42 ml cells + 0.58 ml media
Ala	0.53 ml cells + 0.47 ml media
Raji	0.12 ml cells + 19.88 ml media

40

45

50

55

Results:

Total and Spontaneous CPM

	<u>E:T</u>	<u>Pos. Code</u>	<u>CPM</u>	<u>% Cytolysis</u>
5				
	Blank	1 BLAN	-5.9	.0
		2 BLAN	-8.6	.0
		3 BLAN	-10.3	.0
		0 MEAN	-8.3	.0
10				
	Max.	4 TOTA	823.0	.0
	Release	5 TOTA	700.9	.0
		6 TOTA	896.9	.0
		0 MEAN	806.9	100.0
15				
	Spont.	7 REFR	104.5	13.0
	Release	8 REFR	97.5	12.1
		9 REFR	109.3	13.5
		0 MEAN	103.8	.0
20				

Standard Ficoll 1x10⁶

25				
	20:1	10 UNKS	276.5	24.6
		11 UNKS	287.2	26.1
		12 UNKS	309.1	29.2
		0 MEAN	291.0	26.6
30				
	10:1	13	226.2	17.4
		14	212.4	15.4
		15	255.8	21.6
		0 MEAN	231.4	18.2
35				
	5:1	16	174.4	10.0
		17	180.4	10.9
		18	162.4	8.3
		0 MEAN	172.4	9.8
40				
	2.5:1	19	134.2	4.3
		20	124.8	3.0
		21	99.6	- .6
		0 MEAN	119.5	2.2
45				
50				
55				

Direct - No Ficoll 1×10^6

5	20:1	22	490.6	55.0
		23	436.5	47.3
		24	423.7	45.5
		0 MEAN	450.3	49.3
10	10:1	25	484.5	54.1
		26	526.7	60.1
		27	495.8	55.8
		0 MEAN	502.4	56.7
15	5:1	28	391.6	40.9
		29	395.3	41.5
		30	432.2	46.7
		0 MEAN	406.4	43.0
20	2.5:1	31	236.5	18.9
		32	298.4	27.7
		33	220.7	16.6
		0 MEAN	251.9	21.1

Direct - No Ficoll 5×10^6

30	20:1	34	409.2	43.4
		35	370.9	38.0
		36	361.2	36.6
		0 MEAN	380.5	39.3
35	10:1	37	352.4	35.4
		38	409.3	43.4
		39	364.4	37.1
		0 MEAN	375.3	38.6
40	5:1	40	313.3	29.8
		41	293.3	27.0
		42	292.5	24.0
		0 MEAN	293.1	26.9
45	2.5:1	43	194.5	12.9
		44	227.2	17.6
		45	194.9	13.0
		0 MEAN	205.5	14.5

Direct - No Ficoll 1×10^7

5	20:1	46	245.6	20.2
		47	193.2	12.7
		48	195.7	13.1
		0 MEAN	211.5	15.3
10	10:1	49	191.3	12.4
		50	188.6	12.1
		51	232.2	18.3
		0 MEAN	204.1	14.3
15	5:1	52	174.1	10.0
		53	166.8	9.0
		54	161.8	8.2
		0 MEAN	167.6	9.1
20	2.5:1	55	136.0	4.6
		56	149.8	6.5
		57	114.9	1.6
		0 MEAN	133.6	4.2

25

 ϕ Ala 5×10^6

30	20:1	58	278.1	24.8
		59	221.1	16.7
		60	226.6	17.5
		0 MEAN	241.9	19.6
35	10:1	61	206.5	14.6
		62	188.9	12.1
		63	195.9	13.1
		0 MEAN	197.1	13.3
40	5:1	64	186.5	11.8
		65	199.8	13.7
		66	173.3	9.9
		0 MEAN	186.5	11.8
45	2.5:1	67	199.2	13.6
		68	175.1	10.1
		69	212.7	15.5
		0 MEAN	195.7	13.1

50

EXAMPLE 4

55 A standard leukapheresis product containing 230 mls human leukocytes prepared from 3600 ml whole blood collected in 520 mls anticoagulant ACD-B was obtained from Biological Specialty Corporation, Lansdale, PA. The following procedures were performed using this product.

- 1) Removed 10 mls of cells

A)

- 1) Took 5 mls of this blood and mixed with 5 mls of PBS
- 2) Underlayered 10 mls of Ficoll
- 3) Centrifuged for 30 min @ 2000 rpm
- 4) Washed and counted
- 5) This was the Standard @ 1.5×10^6 cells/ml in 10 ml flask

Standard

Viable = 49

Non-viable = 0

10 % Viable = 100%

Cells/ml = 9.8×10^5 /mlTotal = 1.96×10^8 :20 ml.Dilution: 1.5 ml cells + 8.5 mls media + 1 μ l IL-2

B)

- 15 1) The second 5 mls was used for direct testing
- 2) A WBC (via Unopette®) and differential were performed:
- 3) WBC
 - 4.5×10^7 ml
 - 2.25×10^8 :5 ml
 - 20 Diff.
 - 72% Lymphocytes
 - 22% Granulocytes
 - 6% Monocytes
- 4) Cells were then put up in culture at 1.5×10^6 /ml and 5×10^6 /ml in 10 ml flasks
 - 25 1.5×10^6 /ml = 0.33 ml/cells + 9.67 ml media + 1 μ l IL-2
 - 5×10^6 /ml = 1.11 ml/cells + 8.89 ml media + 1 μ l IL-2
 - Cells were incubated 4 days; chromium release assay was run.

Cell Count:

		Viable	Non-Viable	% Viable	Cells/ml
30	Standard 1.5×10^6	40	3	93	8×10^6
	Direct @ 1.5×10^6	40	2	94	8×10^6
35	Direct @ 5×10^6	24	1	96	4.8×10^6
	Raji	40	6	87	8×10^6

3) LAK Assay

40 E:T ratio 40:1, 20:1, 10:1, 5:1, 2.5:1 1.25:1

Cells were diluted to 4×10^5 Raji's were diluted to 1×10^5

Dilutions

45 Std. 0.5 ml cells + 0.5 ml media

Direct 1.5×10^6 0.5 ml cells + 0.5 ml mediaDirect 5×10^6 0.83 ml cells + 0.17 ml media

Raji 0.25 ml cells + 19.75 ml media

50

55

Results:

Total and Spontaneous CPM			
<u>E:T</u>	<u>Pos. Code</u>	<u>CPM</u>	<u>% Cytolysis</u>
Blank	1 BLAN	-11.0	.0
	2 BLAN	-10.6	.0
	3 BLAN	-10.6	.0
	0 MEAN	-10.7	.0
Max. Release	4 TOTA	576.9	.0
	5 TOTA	564.6	.0
	6 TOTA	584.6	.0
	0 MEAN	575.4	100.0
Spont. Release	7 REFR	139.7	24.3
	8 REFR	140.6	24.4
	9 REFR	152.3	26.5
	0 MEAN	144.2	.0
Standard 1.5x10 ⁶			
40:1	10 UNKS	540.1	91.8
	11 UNKS	547.0	93.4
	12 UNKS	518.3	86.8
	0 MEAN	535.2	90.7
20:1	13	459.5	73.1
	14	503.2	83.3
	15	458.0	72.8
	0 MEAN	473.6	76.4
10:1	16	404.5	60.4
	17	410.8	61.8
	18	439.1	68.4
	0 MEAN	418.1	63.5
5:1	19	320.5	40.9
	20	279.1	31.3
	21	275.0	30.3
	0 MEAN	291.5	34.2
2.5:1	22	222.1	18.1
	23	253.0	25.2
	24	232.1	20.4
	0 MEAN	235.7	21.2
1.25:1	25	175.7	7.3
	26	205.6	14.2
	27	207.9	14.8
	0 MEAN	196.4	12.1

Direct - No Ficoll 1.5×10^6

5	80:1	28	415.5	62.9
		29	393.0	57.7
		30	461.8	73.7
		0 MEAN	423.4	64.8
10	40:1	31	417.7	63.4
		32	404.2	60.3
		33	405.6	60.6
		0 MEAN	409.2	61.5
15	20:1	34	381.1	55.0
		35	407.0	60.9
		36	447.9	70.4
		0 MEAN	412.0	62.1
20	10:1	37	430.5	66.4
		38	407.4	61.0
		39	442.6	69.2
		0 MEAN	426.8	65.5
25	5:1	40	337.5	44.8
		41	357.4	49.5
		42	358.0	49.6
		0 MEAN	351.0	48.0
30	2.5:1	43	240.8	22.4
		44	270.9	29.4
		45	262.1	27.3
		0 MEAN	257.9	26.4
35	1.25:1	46	186.9	9.9
		47	194.1	11.6
		48	182.7	8.9
		0 MEAN	187.9	10.1

Direct - No Ficoll 5×10^6

45	40:1	49	391.0	57.2
		50	392.1	57.5
		51	399.6	59.2
		0 MEAN	394.2	58.0
50	20:1	52	386.3	56.2
		53	381.1	54.9
		54	377.0	54.0
		0 MEAN	381.5	55.0

	10:1	55	346.0	46.8
		56	346.0	46.8
		57	325.5	42.0
5		0 MEAN	339.2	45.2
	5:1	58	296.5	35.3
		59	236.6	21.4
		60	239.0	22.0
10		0 MEAN	257.4	26.2
	2.5:1	61	190.6	10.8
		62	172.8	6.6
		63	190.5	10.7
15		0 MEAN	184.6	9.4
	1.25:1	64	151.7	1.7
		65	172.5	6.6
		66	203.8	13.8
20		0 MEAN	176.0	7.4

EXAMPLE 5

25

Materials:

Buffy coat - 52 mls of blood
 Cell count -
 4.3x10⁷ cells/ml (total WBC)
 1.8x10⁷ neutrophils/ml (42%)
 (est.) 1-2x10⁷ lymphocytes/ml (20-50%)
 (est.) 5x10⁹ RBC/ml
 35 (est.) 0.5-1x10⁷ monocytes/ml
 Ficoll-Hypaque (Ficoll)
 CCM - 5% FCS - RPMI

40 Procedures:1) No Ficoll

45

- a) To 10.5 ml of Buffy coat add 215 ml CCM. Cell Count 2x10⁶ cells/ml (total WBC)
 b) Add 10 μ l/ml of IL-2
 c) Place 112 ml of culture mix in flask
 d) Place 112 ml of culture mix in bag
 50 e) Incubate at 37°C for 20 days.
 f) Sample at 3, 6, 12, 17 and 20 days for cell count and ⁵¹Cr Release (LAK) assay.

2) Ficoll

55

- a) Put 42 mls Buffy Coat in 50 ml centrifuge tube
 b) Centrifuge at 800 g for 10 minutes

c) Discard supernatant, recover mononuclear WBC layer (Lymphocytes and monocytes) floating on Ficoll layer, wash 3X. 300×10^6 total mononuclear cells isolated.

d) Add CCM to provide mononuclear cell concentration of $2 \times 10^6/\text{ml}$

e) Place 75 ml in flask

f) Place 75 ml in bag

g) Incubate at 37°C for 20 days

h) Sample at 3, 6, 12, 17 and 20 days for cell count and ^{51}Cr release (LAK) assay.

Summary of Cell Counts ($\times 10^6/\text{ml}$)

Days Culture	<u>Flasks</u>		<u>Bags</u>	
	<u>Ficoll</u>	<u>No Ficoll</u>	<u>Ficoll</u>	<u>No Ficoll</u>
0	2×10^6	2×10^6	2×10^6	2×10^6
3	1.5×10^6	0.9×10^6	2.1×10^6	$.7 \times 10^6$
6	2×10^6	0.4×10^6	2×10^6	$.8 \times 10^6$
12	2.5×10^6	1.4×10^6	2.7×10^6	2×10^6
17	1.6×10^6	1.2×10^6	1.9×10^6	1.1×10^6
20	1.1×10^6	0.6×10^6	2.4×10^6	1.2×10^6

Summary of LAK Activity

3 LU₃₀

Days Culture	<u>Flasks</u>		<u>Bags</u>	
	<u>Ficoll</u>	<u>No Ficoll</u>	<u>Ficoll</u>	<u>No Ficoll</u>
3	10	40	5	20
6	5	100	5	14
12	2.5	1	2	<1
17	7	<1	7	1
20	7	2	2.5	1

Claims

1. In the method of producing LAK cells *in vitro* which comprises removing RBC's and plasma from whole blood to reduce a lymphocyte-containing WBC-rich fraction and incubating the WBC-rich fraction in culture medium with IL-2, the improvement which comprises removing RBC's and plasma and using the WBC-rich fraction without an intermediate separation of lymphocytes on a ficoll gradient.

2. Method of Claim 1 wherein the RBC's are removed by leukapheresis and the volume percent of RBC's in the WBC-rich fraction is in the range of about 1-20.

3. Method of Claim 2 wherein the RBC/WBC ratio in the WBC-rich fraction is in the range of about 0.2-250.

4. Method of Claim 1 wherein the RBC's are removed by elutriation leukapheresis and the volume percent of RBC's in the WBC-rich fraction is in the range of about 1-6.

5. Method of Claim 4 wherein the RBC/WBC ratio in the WBC-rich fraction is in the range of about 0.2-50.

6. Method of Claim 5 wherein the WBC differential is about 80-85% lymphocytes, about 10-20% monocytes, and about 1-5% granulocytes.

7. Method of Claim 6 wherein the volume percent RBC's in the WBC-rich fraction is in the range of about 2-4, the RBC/WBC ratio in the WBC-rich fraction is in the range of about 0.5-25.

5 8. Method of Claim 1 wherein the monocytes are depleted by treatment with phenyl alanine methyl ester before incubation of the WBC-rich fraction.

9. Method of Claim 2 wherein the WBC-rich fraction is washed with salt solution prior to incubation to inhibit clotting.

10 10. In the method of generating LAK cells by incubating a lymphocyte-containing WBC-rich fraction in culture medium with IL-2, the improvement which comprises using a lymphocyte-containing WBC-rich fraction having a RBC/WBC ratio by number in the range of about 0.2 to 300 and RBC volume percent of about 1-50.

11. Method of Claim 10 wherein the RBC/WBC ratio is in the range of about 0.2-250 and the RBC volume percent is in the range of about 1-20 in the WBC-rich fraction.

15 12. Method of Claim 11 wherein the RBC/WBC ratio is in the range of about 0.2-50 and the RBC volume percent is in the range of about 1-6 in the WBC-rich fraction.

13. Method of Claim 12 wherein the differential of WBC-rich fraction is about 1-5% granulocytes, 0-20% monocytes and greater than about 80% lymphocytes.

20 14. Method of Claim 13 wherein the RBC/WBC ratio is in the range of about 0.5-25 and the RBC volume content is about 2-4 in the WBC-rich fraction.

15 15. In the method of treatment of a cancer patient by adoptive immunotherapy which comprises removing peripheral blood from the patient, separating a lymphocyte-containing WBC-rich fraction from the blood, incubating the lymphocyte-containing WBC-rich fraction with interleukin-2 to produce lymphokine-activated killer cells, and reinjecting the activated cells into the patient, the improvement which comprises separating the lymphocyte-containing WBC-rich fraction without use of a ficoll gradient, whereby the volume percent of RBC's in the WBC-rich fraction is in the range of about 1-20.

16. Method of Claim 15 wherein the lymphocyte-containing WBC-rich fraction is separated by elutriation leukapheresis whereby the volume percent of RBC's in the WBC fraction is in the range of about 1-6.

30

35

40

45

50

55



European Patent
Office

EUROPEAN SEARCH REPORT

Application Number

EP 88 10 6566

DOCUMENTS CONSIDERED TO BE RELEVANT			
Category	Citation of document with indication, where appropriate, of relevant passages	Relevant to claim	CLASSIFICATION OF THE APPLICATION (Int. Cl. 4)
Y,D	US-A-4 464 167 (D.W. SCHOENDORFER et al.) * Whole document *	1-9,15-16	C 12 N 5/00 A 61 K 35/14
Y	PROCEEDINGS OF THE NATIONAL ACADEMIE OF SCIENCES, USA, vol. 82, April 1985, pags 2468-2472; D.L. THIELE et al.: "Regulation of cellular function by products of lysosomal enzyme activity: elimination of human natural killer cells by a dipeptide methyl ester generated from L-leucine methyl ester by monocytes or polymorphonuclear leukocytes" * "Discussion" *	1-9,15,16	
Y	FEDERATION PROCEEDINGS, vol. 44, no. 5, 8th March 1985, 69th annual meeting, Anaheim, California, 21th-26th April 1985, page 1688, abstract no. 7469; T. MEINEKE et al.: "Monocyte modulation of IL-2 induction of lymphokine activated killer cells" * Abstract *	1-9,15,16	
Y	CANCER, vol. 55, 1985, pages 1327-1333; A.A. RAYNER et al.: "Lymphokine-activated killer (LAK) cells" * Abstract *	1-9,15-16	
The present search report has been drawn up for all claims			
Place of search THE HAGUE		Date of completion of the search 02-09-1988	Examiner ALVAREZ Y ALVAREZ C.
CATEGORY OF CITED DOCUMENTS			
X : particularly relevant if taken alone Y : particularly relevant if combined with another document of the same category A : technological background O : non-written disclosure P : intermediate document			
T : theory or principle underlying the invention E : earlier patent document, but published on, or after the filing date D : document cited in the application L : document cited for other reasons ----- & : member of the same patent family, corresponding document			



DOCUMENTS CONSIDERED TO BE RELEVANT			
Category	Citation of document with indication, where appropriate, of relevant passages	Relevant to claim	CLASSIFICATION OF THE APPLICATION (Int. Cl. 4)
Y, D	NEW ENGLAND JOURNAL OF MEDICINE, vol. 313, no. 23, pages 1485-1492, 5th December 1985; S.A. ROSENBERG et al.: "Observations on the systemic administration of autologous lymphokine-activated killer cells and recombinant interleukin-2 to patients with metastatic cancer" * Whole article *	1-9, 15, 16	
			TECHNICAL FIELDS SEARCHED (Int. Cl. 4)
The present search report has been drawn up for all claims			
Place of search THE HAGUE		Date of completion of the search 02-09-1988	Examiner ALVAREZ Y ALVAREZ C.
CATEGORY OF CITED DOCUMENTS X : particularly relevant if taken alone Y : particularly relevant if combined with another document of the same category A : technological background O : non-written disclosure P : intermediate document T : theory or principle underlying the invention E : earlier patent document, but published on, or after the filing date D : document cited in the application L : document cited for other reasons & : member of the same patent family, corresponding document			